2X PCR Solution

PrimeSTAR® HS (Premix)

Code No. R040A
Size: 500 μ l x 5
(for 100 PCR reactions)

Shipping at -20°C
Store at -20°C

Storage: Excessive freezing and thawing may decrease enzyme activity.

Lot No.

Expiration Date:

For a detailed protocol, refer to the product manual supplied with this product.

Description:

PrimeSTAR HS (Premix) is a 2X PCR enzyme mixture composed of PrimeSTAR HS DNA Polymerase, a high-fidelity DNA polymerase, reaction buffer, and dNTP mixture. The premix formulation allows for quick preparation of reactions and reduces the risk of contamination. This product is also useful for high-throughput applications.

PrimeSTAR HS DNA Polymerase is a high-fidelity enzyme that has a very strong 3'-5' exonuclease activity. It also has a higher amplification efficiency than standard *Taq* DNA Polymerase.

Composition:

PrimeSTAR HS DNA Polymerase*: 1.25 units/25 μ l
dNTP Mixture: 2X conc.; ea. 0.4 mM
PrimeSTAR Buffer: 2X conc.; including 2 mM Mg^{2+}

* : Specification of PrimeSTAR HS DNA Polymerase (Cat. #R010A)

Unit definition: One unit is the amount of enzyme that will incorporate 10 nmol of dNTPs into acid-insoluble products in 30 minutes at 74°C with activated salmon sperm DNA as the template-primer.

Reaction mixture for unit definition:

100 mM	Tris-HCl (pH8.3 at 37°C)
10 mM	KCl
2 mM	MgCl_2
6 mM	$(\text{NH}_4)_2\text{SO}_4$
0.1%	Triton X-100
200 μ M	each dATP, dGTP, dCTP
100 μ M	$[\text{^3H}]$ TTP
0.001%	BSA
0.4 mg/ml	activated salmon sperm DNA

Purity: Nicking, endonuclease, and exonuclease activity were not detected after incubation of 0.6 μ g of supercoiled pBR322 DNA, or 0.6 μ g of λ -Hind III digest with 10 units of this enzyme for 1 hour at 74°C .

Applications:

High fidelity DNA amplification by hot start PCR

PCR product:

A significant percentage of PCR products obtained using PrimeSTAR HS DNA polymerase will have blunt-ends. Thus, the obtained PCR products can be directly cloned into blunt-end vectors. If necessary, phosphorylate PCR products before cloning.

PCR test:

Good enzyme performance was confirmed by PCR amplification using both λ DNA (amplified fragments: 8, 10, 12 kb) and human genomic DNA (amplified fragment: 0.5, 1, 2, 4, 6, 8 kb) as templates.

General reaction mixture for PCR (total 50 μ l):

PrimeSTAR HS (Premix)	25 μ l
Primer 1	0.2 - 0.3 μ M (final conc.)
Primer 2	0.2 - 0.3 μ M (final conc.)
Template	< 200 ng
Sterile distilled water	to 50 μ l

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Notice: Because this enzyme has an extremely high priming efficiency, a short annealing time is strongly recommended. Please refer to the user manual supplied with this product for detailed cycling recommendations.

NOTICE TO PURCHASER: LIMITED LICENSE

[P1] PCR Notice

Use of this product is covered by one or more of the following US patents and corresponding patent claims outside the US: 5,789,224, 5,618,711, 6,127,155 and claims outside the US corresponding to expired US Patent No. 5,079,352. The purchase of this product includes a limited, non-transferable immunity from suit under the foregoing patent claims for using only this amount of product for the purchaser's own internal research. No right under any other patent claim, no right to perform any patented method, and no right to perform commercial services of any kind, including without limitation reporting the results of purchaser's activities for a fee or other commercial consideration, is conveyed expressly, by implication, or by estoppel. This product is for research use only. Diagnostic uses under Roche patents require a separate license from Roche. Further information on purchasing licenses may be obtained by contacting the Director of Licensing, Applied Biosystems, 850 Lincoln Centre Drive, Foster City, California 94404, USA.

[M54] PrimeSTAR HS DNA Polymerase

This product is covered by the claims of U.S. Patent No. 7,704,713 and its foreign counterparts.

PrimeSTAR is a registered trademark of TAKARA BIO INC.

Note

This product is for research use only. It is not intended for use in therapeutic or diagnostic procedures for humans or animals. Also, do not use this product as food, cosmetic, or household item, etc. Takara products may not be resold or transferred, modified for resale or transfer, or used to manufacture commercial products without written approval from TAKARA BIO INC. If you require licenses for other use, please contact us by phone at +81 77 543 7247 or from our website at www.takara-bio.com. Your use of this product is also subject to compliance with any applicable licensing requirements described on the product web page. It is your responsibility to review, understand and adhere to any restrictions imposed by such statements. All trademarks are the property of their respective owners. Certain trademarks may not be registered in all jurisdictions.

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